

EDGEBIOCERAMIC™ SEALER

“EdgeBioCeramic Sealer is the best sealer on the market.”

Dr. Charles J. Goodis

Endodontist, Owner, CEO EdgeEndo



EDGEBIOCERAMIC™ SEALER features a unique resin-free formula making it highly biocompatible and hydrophilic. Unlike conventional sealers, Edge BioCeramic Sealer does not shrink and therefore it is not necessary to compact the gutta percha. It can't be any easier...Simply coat the canal walls with EdgeBioCeramic Sealer and place the gutta percha point to length.

OBTURATION HAS BECOME EASIER WITH **EDGEBIOCERAMIC™**.

- Biocompatible
- Premixed
- Hydrophilic
- Highly Antibacterial
- Radiopaque
- Re-treatable
- Zero Shrinkage
- Cost Effective



EDGEENDO®


Twice as Good. Half the Cost.

EdgeEndo.com



1-855-985-3636


Antibacterial (pH) Comparison (Day 1)



HIGHLY ANTIBACTERIAL:

EDGEBIOCERAMIC™ is alkaline (+10pH) making it antibacterial. Compared to conventional sealers EdgeBioceramic Sealer has superior antibacterial properties.*


ISO 6786/2001 Flow Comparison (mm)



FLOW:

EDGEBIOCERAMIC™ exhibits an ideal flow rate which is consistent with ISO 6786/2001.* The small particle size and hydrophilic nature of Edge Bioceramic Sealer allow it to easily flow into canal irregularities.


Calcium Ion Release (Week 1)



BIOCOMPATIBILITY:

EDGEBIOCERAMIC™ actively releases calcium ions.* The unique non-resin bioceramic formula makes Edge Bioceramic highly biocompatible.

Push-Out Bond Strength (Mpa)



BONDING:

EDGEBIOCERAMIC™ features a high push-out bond strength* and forms a tight seal within the root canal thereby reducing the likelihood of bacterial recolonization.

*Studies available upon request