

# EdgeEndo Tip Power Level Recommendations



TIP	Any NEWTRON Device	ProUltra	ProUltra Booster	Ultrawave or Ultrawave XS	Forza	Mini Endo
Edge tip #1	7 (yellow)	7	8	7	E8	40%
Edge tip #2	10 (Yellow)	10	10	10	E10	50%
Edge tip #3	7 (Yellow)	7	8	7	E8	40%
Edge tip #4	7 (Yellow)	7	8	7	E8	40%
Edge tip #5	7 (Yellow)	7	8	7	E8	40%
Edge tip #6	7 (Yellow)	7	8	7	E8	40%
Edge tip #7	7 (Yellow)	7	8	7	E8	40%
Edge tip #8	7 (Yellow)	7	8	7	E8	40%
Edge surgical tip #1	7 (Yellow)	7	8	7	E8	40%
Edge surgical tip #2	7 (Yellow)	7	8	7	E8	40%
Edge surgical tip #3	7 (Yellow)	7	8	7	E8	40%
Edge surgical tip #4	7 (Yellow)	7	8	7	E8	40%
Edge surgical tip #5	7 (Yellow)	7	8	7	E8	40%
Edge surgical tip #6	7 (Yellow)	7	8	7	E8	40%

