



EDGETIP™

EDGESURGICAL TIP™

Tip Compatibility Chart

Tips are compatible	Tips are not compatible
NSK: Varios (any version)	Any EMS or EMS/Hu-Friedy co-branded device
Sybron: Mini Endo	Den Mat: Pro Select Platinum
Obtura Spartan: Wave and MTS	Coltene: Biosonic Suvi
Hu-Friedy: Symmetry IQ	Any Magnetostrictive Device (Cavitron)
Vista: Pilot, Co-Pilot	Kavo: (not ultrasonic)
Bonart: P6 Piezo	Amdent USA: DBI
Woodpecker: Multiple devices (with S handpiece)	Titan: (not ultrasonic)
Dentsply Sirona (Tulsa): Pro Ultra or Pro Ultra Booster	Any device with "E" style threading (E stands for EMS, type)
Vector: Piezo device with S handpiece	Parkell Piezo Tips They are EMS Style threads
A-dec/SP Newtron	Any SIRONA SiroSonic device
Any previous generation ACTEON/SATELEC device: P5 Booster, P5 Newtron (LED), P5 Newtron XS (LED), P Max XS, P Max	
Any device with "S" style threading (S stands for SATELEC=NEWTRON)	
Brasseler Forza V3	